



## Cosyfit Design Pack

The following are instructions on how to produce a design for your Cosyfit.

Please provide contact details, & information about how you intend to use the Cosyfit platform, then complete the Cosyfit Options sheet, to give us details of colours, number of openings & whether a Hi-Lo bed base or tube holes are needed. Our most popular Cosyfit is single bed size with sides 1.47 metres high. We also offer a customised Cosyfit.

Please check that the room where the Cosyfit will be sited has the appropriate clearance around the bed. If you have any questions when completing the room plan, please contact us as we are happy to assist.

Please post, fax or email the completed pack to us at the address below

Safespaces (Cornholme) Ltd  
361 Burnley Road  
Todmorden  
Lancs  
OL14 7DH  
Tel/fax +44 (0)1706 816274  
Email [enquiries@safespaces.co.uk](mailto:enquiries@safespaces.co.uk)



## Contact Details

Details of person requesting quote	
Name	Job title/Role
Address	
Phone	Email
Mobile No	Fax

Details of person using Cosyfit (add separate form if more than one user)	
Parent/Carer's name	Email
Phone	Mobile No
User's Name	Age
Address (location of Safespace)	
Has client had an assessment of need? Yes/No	
Medical condition/s (if any)	
Does client have diagnosis Yes/No State diagnosis	
Mobility  Does user need hoisting? Yes/No  If yes ceiling or mobile hoist?  Can the user pull to stand unassisted?  Can the user walk unaided?	
Behaviours	
Sleep patterns/issues	
Has funding been identified?	

Signature:	Date:
------------	-------


## Cosyfit Design Options

<b>COSYFIT DIMENSIONS</b>	<b>Please tick</b>
Single bed size 91cm width x 183cm length ( suitable for Hi-Lo base)	
Single bed size 91cm width x 198cm length ( suitable for Hi-Lo base)	
Double bed size 122cm width x 183cm length	
Double bed size 122cm width x 198cm length	
King size bed size 183cm width x 183cm length	
King size bed size 183cm width x 198cm length	

<b>HEIGHT OF COSYFIT SIDES.</b> Please note that the mattress is 15cm deep. Please incorporate this into your calculation.	
122cm (can also be used with Hi-lo bases.)	
147cm (can also be used with Hi-lo bases.)	
61cm (for profiling beds)	
76cm (for profiling beds)	

<b>ZIPPED OPENING WITH TRANSPARENT WINDOW (See image on page 6)</b>	
Front zipped opening only	
Both front & rear zipped opening	

<b>SPLIT BACK PANEL</b> We can provide a panel which folds in half if access to your property or to the room where the Cosyfit will be situated is difficult. Please let us know if you think this may be a problem.	
---	--

<b>STANDARD MATTRESS</b> Our standard mattress is covered in strong waterproof PVC and designed to fit snugly into the Cosyfit.		
--	--	--

<b>MEMORY FOAM MATTRESS</b> Our memory foam mattress has a Lycra cover. It is designed to fit snugly into the Safespace and is water resistant. It is not suitable for users who bite objects.	
---	--

<b>TUBE HOLES</b> If these are required, e.g. for oxygen or feeding etc, please state diameter & location of tube hole on diagram on page 6	
--	--

<b>BED BASE OPTIONS</b>	
Choose from a range of bed base heights, & Hi-Lo options below:	
<b>Floor.</b> The Cosyfit sits directly on the floor allowing the user to manoeuvre him/herself in & out	
<b>Platform</b> This is for someone who uses a mobile hoist. The wheels of the hoist can be slid under the platform. The clearance under the platform is 15cm. Please ask if you need a different level of clearance.	
<b>Hi-Lo Trendelenburg</b> This is for someone who needs a hi-lo base to allow for personal care ( with forward & reverse tilt)	
<b>Hi-Lo Profiling</b> bed base (raises knees & head) NB The profiling base is generally ordered with sides up to 76cm.	

<b>MANUAL PROFILING SET</b> For people who need to be supported but who are too mobile to be safe in a profiling bed. Foam wedge for 2 position back support & cylinder roll to act as a knee break	
---	--



When choosing colours, please consider how calm or lively an environment you wish to create for your Cosyfit. The colours below may not be an exact representation, so please check our website. Please choose colours for your Cosyfit, & mattress from the selection below. If there is a particular colour you require please contact us to ask.

COSYFIT COLOUR: \_\_\_\_\_

MATTRESS COLOUR: \_\_\_\_\_



**COSYFIT OPENING PANEL – TUBE HOLES POSITIONS (external)**

As standard the opening panel opens from left to right. If you need the panel to open right to left please tell us.


**Custom Options:** If the design options shown here don't meet your needs. Please call us to discuss your requirements

**HEAD/FOOT**

**END (PLEASE CIRCLE)**

**HEAD/FOOT END**

**(PLEASE CIRCLE)**



**Options for tube holes:** If a Cosyfit **opening back** option has been chosen you can have tube holes in the back. Tube holes will be made 2.5cm diameter unless you state otherwise below. Please tick:

No of tube holes needed:	Head end	Foot end	Diameter needed in cm
One			
Two			
Three			
Four			